



**FOSHAN YEEKONG MEDICAL TECHNOLOGY CO.,LTD.**

---

Add:Luocun Lutang Rd.Shishan,Nanhai,Foshan,PR.China

Tel:+86-757-81260379      86-757-81260787


Email:sales1@yeekongmd.com      export@yeekongmd.com

Website: [www.yeekongmd.com](http://www.yeekongmd.com)

**2023 VERSION**

# Suit Quest Quality

Foshan Yeekong Medical Technology Co., Ltd. is an innovative company specializing in medical rehabilitation device. The company has a professional team and strong production capacity for rehabilitation medical device, the products are manufactured and controlled in accordance with ISO13485 quality regulation. The company will focus on the development of basic and intelligent medical rehabilitation device, including sports rehabilitation equipment, walking aids, smart wheelchairs, and rehabilitation robots. The company will provide customers with innovative and cost-effective products.



# CONTENTS

**03-13**

---

STEEL WHEELCHAIR

**14-23**

---

ALUMINUM WHEELCHAIR

**24-25**

---

POWER WHEELCHAIR

**26-35**

---

COMMUNE

**36-40**

---

BATH BENCH

**41-56**

---

WALKING AIDS

**57**

---

HOSPITAL BED

**58-60**

---

REHABILITATION AID

▶ STEEL WHEELCHAIR



• EK875-46

DIM	94x23x89cm	Max. User Mass( kg)	100
N.W( kg)	16	Main Parameters (cm)	
G.W( kg)	18		



▶ STEEL WHEELCHAIR



• EK909B

DIM	82x28x95cm	Max.User Mass(kg)	100
N.W( kg)	18.2	Main Parameters (cm)	
G.W( kg)	20.5		

▶ STEEL WHEELCHAIR



• EK951B

DIM	83x33x93cm	Max.User Mass(kg)	100
N.W(kg)	21	Main Parameters (cm)	
G.W(kg)	24.1		

▶ STEEL WHEELCHAIR



• EK903AF26

DIM	82x33.5x77cm	Max.User Mass(kg)	100
N.W( kg)	18.2	Main Parameters (cm)	
G.W( kg)	20.5		



▶ STEEL WHEELCHAIR



• EK903ACJF26

DIM	80x33.5x77cm	Max.User Mass(kg)	100
N.W( kg)	18.6	Main Parameters (cm)	
G.W( kg)	20.5		



▶ STEEL WHEELCHAIR



• EK903AJF26

DIM	80x33.5x77cm	Max.User Mass(kg)	100
N.W( kg)	17	Main Parameters (cm)	
G.W( kg)	18.9		

▶ STEEL WHEELCHAIR



• EK809AKD

DIM	80x23x80cm	Max.User Mass(kg)	100
N.W( kg)	16	Main Parameters (cm)	
G.W( kg)	17.7		

▶ STEEL WHEELCHAIR



• EK809F26

DIM	94.5x27.5x89cm	Max.User Mass(kg)	100
N.W( kg)	16.5	Main Parameters (cm)	
G.W( kg)	18.4		



▶ STEEL WHEELCHAIR



• EK903AJF26

DIM	80x33.5x77cm	Max.User Mass(kg)	100
N.W( kg)	17	Main Parameters (cm)	
G.W( kg)	18.9		



## ▶ STEEL WHEELCHAIR



### • EK902C

DIM	80x28x89cm
N.W( kg)	19.3
G.W( kg)	22
Max. User Mass( kg)	100
Main Parameters (cm)	



### • EK869X

DIM	103x25x89cm
N.W( kg)	16.2
G.W( kg)	19
Max. User Mass( kg)	100
Main Parameters (cm)	



### • EK868

DIM	93x28x89cm
N.W( kg)	15
G.W( kg)	17
Max. User Mass( kg)	100
Main Parameters (cm)	



### • EK802-35

DIM	88x23x89cm
N.W( kg)	14.8
G.W( kg)	17
Max. User Mass( kg)	100
Main Parameters (cm)	

## ▶ STEEL WHEELCHAIR



### • EK863ABJ

DIM	76x30x76cm
N.W( kg)	13.1
G.W( kg)	14.4
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK975-51

DIM	80x28x93cm
N.W( kg)	18.7
G.W( kg)	21.3
Max.User Mass( kg)	125
Main Parameters (cm)	



### • EK809J

DIM	94x23x89cm
N.W( kg)	16.7
G.W( kg)	19
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK809Y

DIM	95x31x94cm
N.W( kg)	19
G.W( kg)	21.5
Max.User Mass( kg)	100
Main Parameters (cm)	

▶ ALUMINUM WHEELCHAIR



• EK957LQ

DIM	84x28x92cm	Max.User Mass(kg)	100
N.W( kg)	17.9	Main Parameters (cm)	
G.W( kg)	20.1		



▶ ALUMINUM WHEELCHAIR



• EK908LAJ-46

DIM	80x28x78cm	Max.User Mass(kg)	100
N.W( kg)	15.8	Main Parameters (cm)	
G.W( kg)	17.6		



▶ ALUMINUM WHEELCHAIR



• EK908LAJ-46

DIM	80x28x78cm	Max.User Mass(kg)	100
N.W( kg)	15.8	Main Parameters (cm)	
G.W( kg)	17.6		

▶ ALUMINUM WHEELCHAIR



• EK957LQ

DIM	83x28x91cm	Max.User Mass(kg)	100
N.W( kg)	17.4	Main Parameters (cm)	
G.W( kg)	19.4		

▶ ALUMINUM WHEELCHAIR



• EK018LPQ

DIM	83x32x91cm	Max.User Mass(kg)	100
N.W( kg)	17.1	Main Parameters (cm)	
G.W( kg)	19.1		



## ▶ ALUMINUM WHEELCHAIR



### • EK907LABH-46

DIM	78x28x78cm
NW( kg)	13.1
GW( kg)	15.5
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK908LJ-46

DIM	80x28x91cm
NW( kg)	13.2
GW( kg)	16.3
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK950LBPQ-46

DIM	80x34x93cm
NW( kg)	16.3
GW( kg)	19.4
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK980LA-35

DIM	85x33x73cm
NW( kg)	13
GW( kg)	15.5
Max.User Mass( kg)	100
Main Parameters (cm)	



## ▶ ALUMINUM WHEELCHAIR



### • EK804LABJP-37

DIM	75x25x66cm
NW( kg)	7.2
GW( kg)	9.4
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK805LABJ-46

DIM	71x25x76cm
NW( kg)	8.9
GW( kg)	11.22
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK805LABJP-46

DIM	71x25x76cm
NW( kg)	9
GW( kg)	11.2
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK863LABJPF3-46

DIM	71x29x73.5cm
NW( kg)	10.4
GW( kg)	12.4
Max.User Mass( kg)	100
Main Parameters (cm)	

## ▶ ALUMINUM WHEELCHAIR



### • EK863LABJPF3-46

DIM	71x28x75cm
N.W( kg)	10.6
G.W( kg)	12.4
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK863LAJPF4-46

DIM	80x30x75cm
N.W( kg)	12
G.W( kg)	14.2
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK864L-46

DIM	94x28x92cm
N.W( kg)	12.4
G.W( kg)	15
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK903LQ-46

DIM	80x28x91cm
N.W( kg)	14
G.W( kg)	15.8
Max.User Mass( kg)	100
Main Parameters (cm)	

## ▶ ALUMINUM WHEELCHAIR



### • EKA01-12

DIM	77x57x56cm
N.W( kg)	9.6
G.W( kg)	11
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK863LABJP16-46

DIM	71x29x73.5cm
N.W( kg)	10.4
G.W( kg)	12.4
Max.User Mass( kg)	75
Main Parameters (cm)	



### • EK863LABJF22

DIM	87x36x79cm
N.W( kg)	17
G.W( kg)	19.4
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK958LBCGPY-35/38

DIM	77x40x95cm
N.W( kg)	21
G.W( kg)	24.3
Max.User Mass( kg)	75
Main Parameters (cm)	



## ▶ ALUMINUM WHEELCHAIR



### • EK954LGC-46

DIM	98x35x89cm
N.W( kg)	18.3
G.W( kg)	21
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK985LAF1

DIM	85x39x77cm
N.W( kg)	19.2
G.W( kg)	21.5
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK723LQ-36

DIM	87x36x80cm
N.W( kg)	11
G.W( kg)	13.6
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK778LQ-36

DIM	90x76x68cm
N.W( kg)	12.7
G.W( kg)	14.7
Max.User Mass( kg)	100
Main Parameters (cm)	

## ▶ POWER WHEELCHAIR



### • EK123

DIM	78x68x69cm
N.W( kg)	63
G.W( kg)	70.8
Max.User Mass( kg)	100
km	20
km/h	1~7
Main Parameters (cm)	



### • EK1108

DIM	85.5x35.5x80cm
N.W( kg)	43
G.W( kg)	47.5
Max.User Mass( kg)	120
km	15
km/h	0~6
Main Parameters (cm)	

▶ POWER WHEELCHAIR



• EK1109

DIM	72x66x68cm
N.W(kg)	49
G.W(kg)	55
Max.User Mass(kg)	120
km	20
km/h	0~6
Main Parameters (cm)	



• EK1203

DIM	65x40x80cm
N.W(kg)	28
G.W(kg)	32
Max.User Mass(kg)	120
km	12-15
km/h	0~6
Main Parameters (cm)	





• EK609LBPQUF6-46

DIM	64x61x60cm	Max.User Mass(kg)	100
N.W( kg)	17.7	Main Parameters (cm)	
G.W( kg)	20		

## ▶ COMMODE



### • EK809LPU

DIM	94x29x89cm
NW( kg)	16.6
G.W( kg)	19.6
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK609U

DIM	92x29x89cm
NW( kg)	19
G.W( kg)	21.3
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK681

DIM	75x28x91cm
NW( kg)	19.6
G.W( kg)	21.4
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK609GCU

DIM	107x33x78cm
NW( kg)	25.6
G.W( kg)	29.3
Max.User Mass( kg)	100
Main Parameters (cm)	

▶ **COMMODE**



• **EK691**

DIM	54x53x53cm
NW( kg)	12
G.W( kg)	14
Max.User Mass( kg)	100
Main Parameters (cm)	



• **EK692**

DIM	56x53x57cm
NW( kg)	10.8
G.W( kg)	12.8
Max.User Mass( kg)	100
Main Parameters (cm)	



• **EK698L**

DIM	49x24x67cm
NW( kg)	8.2
G.W( kg)	10
Max.User Mass( kg)	100
Main Parameters (cm)	



• **EK7964L-W5**

DIM	110x35x58cm
NW( kg)	7.1
G.W( kg)	8.6
Max.User Mass( kg)	100
Main Parameters (cm)	



## ▶ COMMODE



### • EK691-2

DIM	58x55x51.5cm
NW( kg)	12
G.W( kg)	14
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK697L

DIM	52x25x70cm
NW( kg)	10.8
G.W( kg)	12.8
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK609BU

DIM	92x29x89cm
NW( kg)	19
G.W( kg)	21.3
Max.User Mass( kg)	100
Main Parameters (cm)	



### • EK609LBU-52

DIM	80x40x93cm
NW( kg)	14
G.W( kg)	16
Max.User Mass( kg)	100
Main Parameters (cm)	

▶ **COMMODE**



• **EK696**

DIM	50x21x76cm
N.W( kg)	7
G.W( kg)	8.2
PCS CN	1
Main Parameters (cm)	



• **EK696L**

DIM	52x21x79cm
N.W( kg)	5.2
G.W( kg)	6.5
PCS CN	1
Main Parameters (cm)	



• **EK813**

DIM	56x40x36cm
N.W( kg)	6
G.W( kg)	7.5
PCS CN	1
Main Parameters (cm)	



• **EK890-2**

DIM	63x32x44cm
N.W( kg)	8
G.W( kg)	8.9
PCS CN	2
Main Parameters (cm)	

▶ **COMMODE**



• **EK893**

DIM	50x19x75cm
N.W( kg)	7.2
G.W( kg)	8.5
PCS CN	1
Main Parameters (cm)	



• **EK894**

DIM	50x21x76cm
N.W( kg)	6
G.W( kg)	7
PCS CN	1
Main Parameters (cm)	



• **EK894L**

DIM	50x23x76cm
N.W( kg)	3.5
G.W( kg)	5
PCS CN	1
Main Parameters (cm)	



• **EK895L**

DIM	60x34x52cm
N.W( kg)	9.2
G.W( kg)	10.6
PCS CN	2
Main Parameters (cm)	



▶ **COMMODE**



• **EK897**

DIM	62X30X42CM
N.W( kg)	9.2
G.W( kg)	10.3
PCS CN	4
Main Parameters (cm)	



• **EK898**

DIM	75x30x40cm
N.W( kg)	9.5
G.W( kg)	10.9
PCS CN	6
Main Parameters (cm)	



• **EK7964L-W3**

DIM	110x35x58cm
N.W( kg)	4.7
G.W( kg)	5.9
PCS CN	1
Main Parameters (cm)	



• **EK8133**

DIM	57x47x52cm
N.W( kg)	7.2
G.W( kg)	8.5
PCS CN	1
Main Parameters (cm)	

▶ **COMMODE**



• **EK8941**

DIM	50x25.5x65cm
N.W( kg)	6.4
G.W( kg)	7.5
PCS CN	1
Main Parameters (cm)	



• **EK8951L**

DIM	45x18x63cm
N.W( kg)	2.6
G.W( kg)	3.6
PCS CN	1
Main Parameters (cm)	



• **EK8101**

DIM	46x22x70cm
N.W( kg)	11.5
G.W( kg)	12.5
PCS CN	2
Main Parameters (cm)	



• **EK811**

DIM	71x52x56cm
N.W( kg)	10
G.W( kg)	11.8
PCS CN	2
Main Parameters (cm)	

▶ **COMMODE**



• **EK8996**

DIM	62x18x80cm
N.W( kg)	7.5
G.W( kg)	8.6
PCS CN	1
Main Parameters (cm)	



• **EK8997-1**

DIM	60x52x93cm
N.W( kg)	15.6
G.W( kg)	17.4
PCS CN	2
Main Parameters (cm)	



• **EK9632L-3**

DIM	58x25x87cm
N.W( kg)	9.5
G.W( kg)	10.7
PCS CN	2
Main Parameters (cm)	



• **EK890-1**

DIM	72x25x39cm
N.W( kg)	7.2
G.W( kg)	8.5
PCS CN	2
Main Parameters (cm)	



▶ **COMMODE**



• **EK899L-1**

DIM	55x20x73cm
N.W( kg)	5.74
G.W( kg)	6.94
PCS CN	1
Main Parameters (cm)	



• **EK899-A**

DIM	79x43x55cm
N.W( kg)	16.7
G.W( kg)	18.7
PCS CN	4
Main Parameters (cm)	



• **EK696L-1**

DIM	52x21x73cm
N.W( kg)	5.48
G.W( kg)	6.7
Max. User Mass( kg)	100
Main Parameters (cm)	



• **EK4002**

DIM	71x52x56cm
N.W( kg)	10
G.W( kg)	11.8
PCS CN	1
Main Parameters (cm)	

## ▶ BATH BENCH



### • EK568S

DIM	45x27x19cm
N.W( kg)	3.2
G.W( kg)	3.7
PCS CN	2
Main Parameters (cm)	



### • EK569S

DIM	50x42x46cm
N.W( kg)	6.4
G.W( kg)	7.6
PCS CN	2
Main Parameters (cm)	



### • EK790

DIM	74x19x52cm
N.W( kg)	6.4
G.W( kg)	7.2
PCS CN	2
Main Parameters (cm)	



### • EK794L

DIM	78x52x23cm
N.W( kg)	8.1
G.W( kg)	9.3
PCS CN	5
Main Parameters (cm)	

## ▶ BATH BENCH



### • EK797L

DIM	60x17x36.5cm
N.W( kg)	3.6
G.W( kg)	4.4
PCS CN	2
Main Parameters (cm)	



### • EK798L

DIM	68x19x39cm
N.W( kg)	5.9
G.W( kg)	7.1
PCS CN	2
Main Parameters (cm)	



### • EK7954

DIM	57.5x39x25cm
N.W( kg)	8.1
G.W( kg)	9
PCS CN	2
Main Parameters (cm)	



### • EK7972L

DIM	35x35x39cm
N.W( kg)	3.1
G.W( kg)	3.9
PCS CN	2
Main Parameters (cm)	



## ▶ BATH BENCH



### • EK7981L

DIM	51x47x48cm
N.W( kg)	3
G.W( kg)	4
PCS CN	1
Main Parameters (cm)	



### • EK7982L

DIM	68x19x39cm
N.W( kg)	6
G.W( kg)	7
PCS CN	2
Main Parameters (cm)	



### • EK7985L

DIM	56x52x45cm
N.W( kg)	7.6
G.W( kg)	8.2
PCS CN	2
Main Parameters (cm)	



### • EK7985LA

DIM	56x52x45cm
N.W( kg)	7.6
G.W( kg)	8.2
PCS CN	2
Main Parameters (cm)	

## ▶ BATH BENCH



### • EK7992LU

DIM	74x51x91.5cm
NW(kg)	11
G.W(kg)	13.5
PCS CN	2
Main Parameters (cm)	



### • EK7985LQ-1

DIM	55x25x36cm
NW(kg)	7.2
G.W(kg)	8.5
PCS CN	2
Main Parameters (cm)	



### • EK7981LQ

DIM	51x47x43cm
NW(kg)	3.2
G.W(kg)	4.2
PCS CN	1
Main Parameters (cm)	




### • EK4001

DIM	35x35x39cm
NW(kg)	3.1
G.W(kg)	3.9
PCS CN	2
Main Parameters (cm)	

## ▶ BATH BENCH




### • EK7931-1

DIM	32x7x71.5cm
N.W( kg)	2.8
G.W( kg)	3.3
PCS CN	1
Main Parameters (cm)	




### • EK791S

DIM	70x56x41.5cm
N.W( kg)	8.5
G.W( kg)	10.6
PCS CN	2
Main Parameters (cm)	




### • EK7986L-B

DIM	62x15x42cm
N.W( kg)	3.1
G.W( kg)	3.9
PCS CN	1
Main Parameters (cm)	



### • EK7972L-A


DIM	35x8x46.5cm
N.W( kg)	1.7
G.W( kg)	2.1
PCS CN	1
Main Parameters (cm)	



## ▶ WALKING AIDS




### • EK965LH-1

DIM	58x22x54cm
N.W( kg)	7.5
G.W( kg)	8.5
PCS CN	1
Main Parameters (cm)	




### • EK965LH-2

DIM	62x24x84cm
N.W( kg)	7.4
G.W( kg)	8.1
PCS CN	1
Main Parameters (cm)	




### • EK965LHF7

DIM	42x28x59cm
N.W( kg)	9.65
G.W( kg)	10.65
PCS CN	1
Main Parameters (cm)	



### • EK965LHF8

DIM	42x68x64cm
N.W( kg)	8.2
G.W( kg)	9.2
PCS CN	1
Main Parameters (cm)	

## ▶ WALKING AIDS



### • EK965LH11

DIM	61x23x85cm
N.W( kg)	7.5
G.W( kg)	9.1
PCS CN	1
Main Parameters (cm)	



### • EK965LH12

DIM	58x22x54cm
N.W( kg)	7.5
G.W( kg)	8.5
PCS CN	1
Main Parameters (cm)	



### • EK965LH-D

DIM	52x18x55cm
N.W( kg)	5.86
G.W( kg)	7
PCS CN	1
Main Parameters (cm)	



### • EK9650LH

DIM	61x23x86cm
N.W( kg)	8.8
G.W( kg)	10.4
PCS CN	1
Main Parameters (cm)	

## ▶ WALKING AIDS



### • EK965LH

DIM	62x24x84cm
N.W( kg)	7.4
G.W( kg)	9
PCS CN	1
Main Parameters (cm)	



### • EK966LH

DIM	62x24x84cm
N.W( kg)	7
G.W( kg)	8.5
PCS CN	1
Main Parameters (cm)	



### • EK965LH4-1

DIM	47.5x33x36.5cm
N.W( kg)	9.6
G.W( kg)	10.2
PCS CN	1
Main Parameters (cm)	



### • EK965LH4-2

DIM	52x22x44cm
N.W( kg)	9.5
G.W( kg)	10.15
PCS CN	1
Main Parameters (cm)	



## ▶ WALKING AIDS



### • EK965LH4-3

DIM	60x33x35cm
N.W( kg)	10.25
G.W( kg)	11
PCS CN	1
Main Parameters (cm)	



### • EK965LH4-4

DIM	56x45x38cm
N.W( kg)	9.75
G.W( kg)	10.95
PCS CN	1
Main Parameters (cm)	



### • EK965LH4-5

DIM	52x21x41cm
N.W( kg)	8.75
G.W( kg)	9.65
PCS CN	1
Main Parameters (cm)	



### • EK965LH3-1

DIM	60x33x35cm
N.W( kg)	11.2
G.W( kg)	12
PCS CN	1
Main Parameters (cm)	

## ▶ WALKING AIDS



### • EK203

DIM	59x32x84cm
N.W( kg)	18
G.W( kg)	20
PCS CN	1
Main Parameters (cm)	



### • EK201

DIM	59x32x84cm
N.W( kg)	20
G.W( kg)	22
PCS CN	1
Main Parameters (cm)	



### • EK9610

DIM	75x15x92cm
N.W( kg)	11.9
G.W( kg)	13.6
PCS CN	1
Main Parameters (cm)	




### • EK9612

DIM	75x22x96cm
N.W( kg)	9.4
G.W( kg)	11.8
PCS CN	1
Main Parameters (cm)	

## ▶ WALKING AIDS




### • EK912L

DIM	58x18x82cm
NW( kg)	5.5
GW( kg)	6.7
PCS CN	2
Main Parameters (cm)	




### • EK912L(XS)-3

DIM	53x19x60cm
NW( kg)	4.3
GW( kg)	5.5
PCS CN	2
Main Parameters (cm)	




### • EK914L

DIM	60x34x89.5cm
NW( kg)	7
GW( kg)	8.7
PCS CN	2
Main Parameters (cm)	



### • EK9122L

DIM	66x18x65cm
NW( kg)	6
GW( kg)	6.9
PCS CN	2
Main Parameters (cm)	

## ▶ WALKING AIDS



### • EK9122L-A

DIM	61x25x80cm
N.W( kg)	8.6
G.W( kg)	9.7
PCS CN	2
Main Parameters (cm)	



### • EK9123L(S)-B

DIM	68x29x59cm
N.W( kg)	6.6
G.W( kg)	7.9
PCS CN	2
Main Parameters (cm)	



### • EK9121L-A

DIM	59x30x52cm
N.W( kg)	5.2
G.W( kg)	6.2
PCS CN	2
Main Parameters (cm)	



### • EK9125L

DIM	62x18x84cm
N.W( kg)	6
G.W( kg)	6.9
PCS CN	2
Main Parameters (cm)	



## ▶ WALKING AIDS



### • EK9632L-1

DIM	66x25x81cm
N.W( kg)	8.2
G.W( kg)	9.5
PCS CN	2
Main Parameters (cm)	



### • EK912L(XS)

DIM	53x18x52cm
N.W( kg)	3.6
G.W( kg)	4.6
PCS CN	2
Main Parameters (cm)	



### • EK915L(XS)-1

DIM	53x19x60cm
N.W( kg)	4.8
G.W( kg)	6
PCS CN	2
Main Parameters (cm)	



### • EK915L(XS)-3

DIM	53x19x65cm
N.W( kg)	6
G.W( kg)	6.9
PCS CN	2
Main Parameters (cm)	

## ▶ WALKING AIDS



### • EK919L

DIM	53x17x81cm
NW( kg)	4.8
GW( kg)	5.8
$\frac{PCS}{CN}$	2
Main Parameters (cm)	



### • EK961L

DIM	57x19x79cm
NW( kg)	6.7
GW( kg)	7.7
$\frac{PCS}{CN}$	2
Main Parameters (cm)	



### • EK963L

DIM	58x18x81cm
NW( kg)	5.2
GW( kg)	7
$\frac{PCS}{CN}$	2
Main Parameters (cm)	



### • EK967L

DIM	57x18x82cm
NW( kg)	5.2
GW( kg)	6.5
$\frac{PCS}{CN}$	2
Main Parameters (cm)	

## ▶ WALKING AIDS



### • EK913L

DIM	53x17x81cm
N.W( kg)	4.8
G.W( kg)	5.7
$\frac{PCS}{CN}$	2
Main Parameters (cm)	



### • EK915L

DIM	53x17x81cm
N.W( kg)	5.1
G.W( kg)	6.2
$\frac{PCS}{CN}$	2
Main Parameters (cm)	



### • EK917L(S)

DIM	76x68x35cm
N.W( kg)	14.9
G.W( kg)	16.9
$\frac{PCS}{CN}$	10
Main Parameters (cm)	

### • EK917L(M)

DIM	86x68x35cm
N.W( kg)	15.2
G.W( kg)	17.5
$\frac{PCS}{CN}$	10
Main Parameters (cm)	


### • EK917L(L)

DIM	96x68x35cm
N.W( kg)	19
G.W( kg)	22.1
$\frac{PCS}{CN}$	10
Main Parameters (cm)	

## ▶ WALKING AIDS




### • EK9145L

DIM	52.5x22x63cm
N.W( kg)	2.9
G.W( kg)	3.9
<u>PCS</u> CN	2
Main Parameters (cm)	




### • EK9637L

DIM	57x19x86cm
N.W( kg)	5.1
G.W( kg)	6.4
<u>PCS</u> CN	2
Main Parameters (cm)	




### • EK9635L

DIM	57x19x86cm
N.W( kg)	5.1
G.W( kg)	6.4
<u>PCS</u> CN	2
Main Parameters (cm)	



### • EK913L-1


DIM	54x25x82cm
N.W( kg)	8.8
G.W( kg)	10
<u>PCS</u> CN	2
Main Parameters (cm)	



## ▶ WALKING AIDS




### • EK915L(XS)

DIM	53x19x60cm
N.W( kg)	4.1
G.W( kg)	5.3
PCS CN	2
Main Parameters (cm)	




### • EK9632L-3

DIM	58x25x86cm
N.W( kg)	7.8
G.W( kg)	9
PCS CN	2
Main Parameters (cm)	




### • EK8143

DIM	47x15x34cm
N.W( kg)	2.8
G.W( kg)	3.6
PCS CN	2
Main Parameters (cm)	



### • EK8511

DIM	39x9x46cm
N.W( kg)	2.5
G.W( kg)	3
PCS CN	2
Main Parameters (cm)	

## ▶ WALKING AIDS



### • EK8512

DIM	53x10x58cm
N.W( kg)	6.1
G.W( kg)	7
PCS CN	2
Main Parameters (cm)	



### • EK9152L

DIM	53x17x81cm
N.W( kg)	5.1
G.W( kg)	6.2
PCS CN	2
Main Parameters (cm)	



### • EK9503L

DIM	50x17x60cm
N.W( kg)	4.5
G.W( kg)	5.4
PCS CN	2
Main Parameters (cm)	



### • EK9126L

DIM	59x18x70cm
N.W( kg)	5
G.W( kg)	6.1
PCS CN	2
Main Parameters (cm)	

# ▶ WALKING AIDS



**EK910L**

DIM	94x15x28cm
Hight(cm)	71-94
N.W( kg)	0.35
G.W( kg)	7.9
PCS-CN	20



**EK920L**

DIM	65x15x26cm
Hight(cm)	71-94
N.W( kg)	0.36
G.W( kg)	7.7
PCS-CN	20



**EK927L**

DIM	62x16x26cm
Hight(cm)	80-90
N.W( kg)	0.37
G.W( kg)	8.3
PCS-CN	20



**EK929L**

DIM	65x15x26cm
Hight(cm)	71-94
N.W( kg)	0.36
G.W( kg)	7.8
PCS-CN	20



**EK930L**

DIM	65x15x27cm
Hight(cm)	71-94
N.W( kg)	0.35
G.W( kg)	7.6
PCS-CN	20



**EK939L**

DIM	63x15x28cm
Hight(cm)	67-90
N.W( kg)	0.34
G.W( kg)	7.3
PCS-CN	20



**EK9206L**

DIM	65x15x26cm
Hight(cm)	71-94
N.W( kg)	0.33
G.W( kg)	7
PCS-CN	20



**EK928L**

DIM	79x29x18cm
Hight(cm)	77-100
N.W( kg)	0.35
G.W( kg)	7.9
PCS-CN	20



**EK938L**

DIM	79x29x18cm
Hight(cm)	77-100
N.W( kg)	0.32
G.W( kg)	7.5
PCS-CN	20



**EK9280L**

DIM	68x35x18cm
Hight(cm)	75-98
N.W( kg)	0.33
G.W( kg)	7.4
PCS-CN	20



**EK9281L**

DIM	68x35x18cm
Hight(cm)	75-98
N.W( kg)	0.33
G.W( kg)	7.5
PCS-CN	20



**EK9300L**

DIM	74x33x14cm
Hight(cm)	75-98
N.W( kg)	0.4
G.W( kg)	8.8
PCS-CN	20

# ▶ WALKING AIDS



**EK923L**

DIM	86x26x29cm
Hight(cm)	86-117
N.W( kg)	0.52
G.W( kg)	11.3
PCS-CN	20



**EK937L**

DIM	100x30x30cm
Hight(cm)	96-119
N.W( kg)	0.53
G.W( kg)	11.5
PCS-CN	20



**EK9332L**

DIM	95x29x31cm
Hight(cm)	95-126
N.W( kg)	0.57
G.W( kg)	12.1
PCS-CN	20



**EK925L(L)**

DIM	140x33x29cm
Hight(cm)	135-155
N.W( kg)	0.88
G.W( kg)	19.7
PCS-CN	20



**EK925L(M)**

DIM	119x33x29cm
Hight(cm)	115-135
N.W( kg)	0.82
G.W( kg)	18.4
PCS-CN	20



**EK925L(S)**

DIM	101x33x29cm
Hight(cm)	95-115
N.W( kg)	0.73
G.W( kg)	16.7
PCS-CN	20



**EK926**

DIM	78x29x36cm
Hight(cm)	74-97
N.W( kg)	0.76
G.W( kg)	8.7
PCS-CN	10



**EK921**

DIM	79x34x44cm
Hight(cm)	72-95
N.W( kg)	0.87
G.W( kg)	10.4
PCS-CN	10



**EK922**

DIM	78x29x43cm
Hight(cm)	74-97
N.W( kg)	0.83
G.W( kg)	9.5
PCS-CN	10



**EK924**

DIM	72x30x34cm
Hight(cm)	65-88
N.W( kg)	0.8
G.W( kg)	9.2
PCS-CN	10



**EK931**

DIM	79x34x44cm
Hight(cm)	72-95
N.W( kg)	0.9
G.W( kg)	10.6
PCS-CN	10



**EK932**

DIM	78x29x43cm
Hight(cm)	74-97
N.W( kg)	0.88
G.W( kg)	9.7
PCS-CN	10



# ▶ WALKING AIDS



**EK934**

DIM	72x30x34cm
Hight(cm)	65-88
N.W( kg)	0.78
G.W( kg)	9.1
PCS-CN	10

**EK945**

DIM	79x34x44cm
Hight(cm)	73-88
N.W( kg)	0.78
G.W( kg)	8.8
PCS-CN	10

**EK946S**

DIM	78x29x43cm
Hight(cm)	74-97
N.W( kg)	0.85
G.W( kg)	9.5
PCS-CN	10

**EK947**

DIM	78x29x36cm
Hight(cm)	72-95
N.W( kg)	0.77
G.W( kg)	8.8
PCS-CN	10

**EK9207LW-1**

DIM	65x15x26cm
Hight(cm)	71-94
N.W( kg)	0.36
G.W( kg)	7.7
PCS-CN	20

**EK933L-4**

DIM	98x33x48cm
Hight(cm)	90-113
N.W( kg)	0.98
G.W( kg)	20.5
PCS-CN	20

**EK941**

DIM	81x36x45cm
Hight(cm)	72-95
N.W( kg)	0.9
G.W( kg)	10.7
PCS-CN	10



**EK942S**

DIM	78x29x43cm
Hight(cm)	74-97
N.W( kg)	0.85
G.W( kg)	9.7
PCS-CN	10

**EK9322**

DIM	78x29x43cm
Hight(cm)	72-97
N.W( kg)	0.81
G.W( kg)	9.3
PCS-CN	10

**EK946A**

DIM	78x29x43cm
Hight(cm)	74-97
N.W( kg)	0.85
G.W( kg)	9.5
PCS-CN	10

**EK921L**

DIM	83x41x46cm
Hight(cm)	74-97
N.W( kg)	0.74
G.W( kg)	9.4
PCS-CN	10

**EK931L**

DIM	83x41x46cm
Hight(cm)	74-97
N.W( kg)	0.74
G.W( kg)	9.4
PCS-CN	10

**EK911L**

DIM	85x26x40cm
Hight(cm)	84
N.W( kg)	0.87
G.W( kg)	9.7
PCS-CN	10


**EK911L-1**

DIM	86x30x41cm
Hight(cm)	84-94
N.W( kg)	1.0
G.W( kg)	11.5
PCS-CN	10



## ▶ HOSPITAL BED



### • EK3230W



### • EK3238W



▶ REHABILITATION AID



• EK5773



• EK302



• EK5777



• EK9601



• EK960



• EK578

L:45cm W:47cm H:78cm



• EK578B

L:45cm W:47cm H:78cm



• EK505

▶ REHABILITATION AID



• EK562

L:77cm W:40cm H:72-113cm



• EK572

L:77cm W:40cm H:72-113cm



• EK5623

L:90cm W:45cm H:72-113cm



• EK5721

L:77cm W:38cm H:72-112cm



• EK5625

L:90cm W:45cm H:72-113cm



• EK531



▶ REHABILITATION AID



• EK8141



• EK661B



• EK662B



• EK663B



• EK664B



• EK665B



• EK115 (01)



• EK115 (02)



• EK115 (03)



• EK667B



• EK666B-100/50



• EK8142



**FOSHAN YEEKONG MEDICAL TECHNOLOGY CO.,LTD.**

---

Add:Luocun Lutang Rd.Shishan,Nanhai,Foshan,PR.China

Tel:+86-757-81260379      86-757-81260787

Email:sales1@yeekongmd.com      export@yeekongmd.com

Website: [www.yeekongmd.com](http://www.yeekongmd.com)